

Scale-up and manufacturing comparability of an Sf RVN baculovirus production process from 10 L to 200 L for AAV

Nene Kalu¹, Eric Lin¹, Noah Miller-Medzon¹, Christina Pelella¹, Allie Van Pelt¹, Aileen Rottinger¹, Harris Shaikh¹, Eric Tompkins¹, Rebecca Bamidele¹, Elena Bianchetti¹, José Manuel Otero¹, Timothy Fenn¹, Vincent Turula¹

¹Lexeo Therapeutics, New York, NY, United States



INTRODUCTION

- Robust, scalable manufacturing processes are critical to supporting late-stage clinical development and eventual commercialization of recombinant adeno-associated virus (rAAV) gene therapies¹
- Traditional transient transfection-based mammalian platforms are often limited by high material costs, low yields, and scalability constraints, particularly for programs requiring high doses or large patient populations¹
- To address these challenges, Lexeo Therapeutics has developed a scalable Sf9 baculovirus-infected insect cell (BIC) platform to produce adeno-associated virus serotype rh.10 (AAVrh.10) vectors
- This study evaluates process scale-up and manufacturing comparability of an Sf9 baculovirus production process for rAAV production, assessing whether a 10 L small-scale model is representative and predictive of a 200 L production-scale stirred-tank bioreactor (STR) process

METHODS

Manufacturing Strategy

- Three process development drug substance (DS) batches (two at 10 L scale and one at 200 L scale) and one 200 L scale Good Manufacturing Practice DS batch were manufactured using conserved critical starting materials and aligned process conditions, where feasible


Upstream Process (Figure 1)

- A rhabdovirus-negative *Spodoptera frugiperda* (Sf-RVN) suspension cell line was cultured in chemically defined, animal-component-free medium
- Two GLP recombinant baculovirus (rBV) constructs were utilized: one containing the *Rep* and *Cap* genes for the AAVrh.10 capsid, and the other containing the gene of interest (GOI)
- The Sf9 cell line was infected with rBV to generate AAVrh.10 and GOI BICs
- For rAAV production, the Sf9 master cell bank thawed and expanded in serial shake flasks, followed by inoculation of the STR and infection at low multiplicity of infection with BICs
- Bioreactor parameters such as agitation and gassing were scaled using power density (P/V) and volumetric gas flow (VVM), respectively

Downstream Process (Figure 1)

- Harvested material underwent chemical lysis and nuclease treatment, followed by depth filtration and tangential flow filtration (TFF)-1
- Purification was performed using affinity (AFF) chromatography and anion exchange (AEX) chromatography
- Final concentration and buffer exchange were completed by TFF-2, followed by sterile filtration to generate bulk drug substance (BDS)

Figure 1: Sf9 baculovirus manufacturing process



Steps indicated in the figure refer to: (1) Sf9-RVN thaw; (2) cell expansion; (3) STR inoculation, infection, and feed addition; (4) harvest and clarification; (5) TFF-1; (6) AFF and (7) AEX chromatography; (8) TFF-2; (9) filtration and BDS fill; and (10) final drug product.

AEX, anion exchange; AFF, affinity; BDS, bulk drug substance; Sf9-RVN, rhabdovirus-negative *Spodoptera frugiperda*; STR, stirred-tank bioreactor; TFF, tangential flow filtration

Comparability Assessment

- Comparability across scales was evaluated using:
 - In-process controls (cell growth, metabolites [data not shown], infection kinetics)
 - Step and overall downstream recoveries
 - Chromatographic performance
 - Critical quality attributes (CQAs), including genome and capsid titer, full/empty capsid distribution, capsid purity, viral protein ratios, and vector genome integrity

RESULTS

Upstream Performance

- Characteristics such as cell growth, viability, and cell diameter were comparable across scales
- Infection kinetics showed consistent asynchronous infection behavior at 10 L and 200 L scales (Figure 2)
- Genome titers (ddPCR) and capsid titers (ELISA) measured in crude viral lysate demonstrated consistent productivity and packaging efficiency during scale-up from 10 L to 200 L bioreactor scales (Figure 3)

Downstream Performance

- Depth filtration, TFF, AFF and AEX chromatography all performed reproducibly at 10 L and 200 L scales
- AFF and AEX chromatograms demonstrated comparable peak shapes and elution profiles across scales
- AEX chromatography achieved consistent enrichment of full capsids, with $\leq 20\%$ empty capsids in BDS (Figure 4)

Yield and Recovery

- Comparable step recoveries were observed across process scales, indicating robust scale-up of downstream AAV unit operations (Figure 5)
- Overall downstream recovery was $\geq 50\%$ at both 10 L and 200 L scales, with higher recoveries observed at 200 L, consistent with improved manufacturing efficiency at larger processing volumes
- Purified yields reached $>1 \times 10^{15}$ vg/L at the 200 L scale, supporting high-dose and large-population clinical supply

Figure 2: Cell culture during AAV production at 10 L and 200 L bioreactor scales




Figure 3: Crude lysate genome and capsid titers at 10 L and 200 L bioreactor scales




Figure 4: Enrichment of full capsids across the downstream purification unit operations*





Figure 5: Step recoveries by ddPCR at 10 L and 200 L bioreactor scales*



PRODUCT QUALITY

- Final DS from 10 L and 200 L scale processes showed comparable CQAs (Table 1):
 - Capsid purity of $\sim 97-98\%$
 - Consistent VP1:VP2:VP3 ratios
 - Comparable full/empty capsid distributions by mass photometry and sedimentation velocity analytical ultracentrifugation
 - Expected vector genome length by long-read sequencing

Table 1. Product quality characteristics from runs at 10 L and 200 L bioreactor scales

	Capsid Purity	VP1:VP2:VP3	Empty Capsids (by SV-AUC)	Genome Sequence Length/Expected
10 L Bioreactor A	98.3%	2.0:1.3:10	12.9%	Pending
10 L Bioreactor B	98.4%	2.3:1.3:10	15.1%	99.9%
200 L PD	98.2%	2.0:1.2:10	19.2%	99.9%
200 L GMP	97.5%	1.9:1.5:10	19.4%	Pending

GMP, good manufacturing practice; PD, process development; SV-AUC, sedimentation velocity analytical ultracentrifugation; VP, viral protein

CONCLUSION

- The successful scale-up of Lexeo Therapeutics' Sf9 baculovirus rAAV production process from 10 L to 200 L scale was demonstrated
- The 10 L small-scale model is representative and predictive of the 200 L production-scale STR process
- Comparable productivity, downstream performance, and product quality were maintained across scales
- This scalable platform supports process development, validation strategy planning, and clinical and commercial manufacturing of AAVrh.10 gene therapies

KEY TAKEAWAYS

- Lexeo Therapeutics' Sf9 baculovirus platform overcomes key limitations of plasmid-based rAAV manufacturing by minimizing the use of critical starting materials
- High purified yields ($>10^{15}$ vg/L), strong overall recoveries ($\geq 50\%$), and favorable empty capsid profiles ($\leq 20\%$) were achieved at the 200 L production scale
- Demonstrated scalability positions this process for robust commercialization and improved patient access to gene therapies

MORE INFORMATION



Learn more about Lexeo



Read our previous publications

REFERENCE

1. Fu Q, et al. *Biotechnol Bioeng.* 2023;120(9):2601-2621.

ACKNOWLEDGMENTS

We acknowledge the technical support and feedback provided by the Lexeo Therapeutics CMC team, the analytical support provided by our CMC Analytical team, and poster support from the Medical Affairs team. We acknowledge technical collaboration and support provided by our colleagues at Millipore and Virovek, and our manufacturing partners at Fujifilm. Editorial support was provided by Sharmin Bovill (BOLDSCIENCE Ltd. UK) and was funded by Lexeo Therapeutics. This poster was developed in accordance with Good Publication Practice (GPP) guidelines.

DISCLOSURES

All authors are employees of Lexeo Therapeutics, Inc., New York, NY, United States.